

Femtosecond Time-Resolved Infrared-Resonant Third-Order Sum-Frequency Spectroscopy

Jizhou Wang^{1,2}, Kai Wang^{1,3*}, Yujie Shen¹, Zehua Han¹, Fu Li¹, Zhe He^{1,4},
Da-wei Wang^{2,5}, Alexei V. Sokolov^{1,6*}, Marlan O. Scully^{1,6}

¹Institute for Quantum Science and Engineering, Texas A&M University, College Station, Texas 77840, USA

²Interdisciplinary Center for Quantum Information and State Key Laboratory of Modern Optical Instrumentation, Zhejiang Province Key Laboratory of Quantum Technology and Device, and Department of Physics, Zhejiang University, Hangzhou, Zhejiang 310027, China

³School of Physics and Astronomy, Sun Yat-sen University, Zhuhai, Guangdong, 519082, China

⁴Department of Medical Engineering, California Institute of Technology, Pasadena, California 91125, USA

⁵Zhejiang Laboratory, Hangzhou, Zhejiang 311121, China

⁶Baylor Research Innovative Center, Baylor University, Waco, Texas 76798, USA

S1- 2D Spectrum Simulation

We simulate the 2D spectrum of the film by considering the strong near-resonant absorption and pre-pulse effects. ITS spectroscopy is sensitive to the IR active modes. For the sample we use, they are at 2853 cm^{-1} and 2926 cm^{-1} separately according to the FTIR spectrum. The full width at half maximum (FWHM) of the 2853 cm^{-1} and 2926 cm^{-1} absorption lines are set as 40 cm^{-1} and 95 cm^{-1} , respectively. In this simulation, we assume that the probe NIR pulse contained a pre-pulse, which is π phase difference from the main pulse. The pre pulse is due to the edge bandpass filter, which results in a sinc function for the pulse. The simulated results are shown in fig. S1.




Figure S1. (a) Simulated 2D spectrum of the time-resolved ITS signal generated from the film. The X-axis is the delay time between the NIR and MIR pulses. The zero probe delay is the time

when the two pulses overlaps. Y-axis is the wavenumber of the visible signal. (b) Signal intensity changes with the time delay at 15740 cm^{-1} . X-axis is the time delay. Y-axis is the intensity.