


## Supplementary Information




**SFig1. Nanoparticle tracking analysis (NTA) corresponding to sample shown in Figure 4.** There is a general broad peak corresponding to particles 30-50 nm in diameter, which is similar to the ~44 nm vesicle size we report with FLOWER; however, there are also two peaks of larger particle size. In contrast to FLOWER, NTA does not have antibody-functionalization for specific exosome detection. Based on our results we expect particles of a larger size found in the Nanosight do not specifically bind to the toroid and do not represent the specific exosomes we are trying to detect.

### Dielectric factor (D) calculation:


Our dielectric factor is defined as:

$$D = \frac{4\pi n_m^2 (n_p^2 - n_m^2)}{(n_p^2 + 2n_m^2)}$$


where  $n_m$  is the index of refraction of the surrounding media, and  $n_p$  is the index refraction of the particle. Using 1.375 as the index of refraction of an exosome, and 1.33 for the index of refraction of the surrounding media, we calculate a  $D \cong 0.5$ .



**SFig2. Resonance wavelength as a function of time for serum from a mouse with no palpable tumor.** Serum from a mouse with no tumor taken from week 5 does not generate a larger signal than serum taken from week 1. This is in contrast to samples taken from mice containing tumors that are analyzed using FLOWER.



**SFig3. Zoom-in of week 4 (Figure 2).** Discrete steps may also be observed upon zooming in on week 4 (Figure 2).



**SFig4. Control performed with an antibody other than anti-CD81.** (a) Resonance wavelength as a function of time for the same serum shown in Fig. 4a flowed over a resonator functionalized with monoclonal anti-tobacco mosaic virus (TMV) instead of anti-CD81 (nine traces). (b) A closer look at a region from (a). No permanent binding steps of exosome-sized particles are observed.