

*Supporting Information*  
*for*

**Fluorogenic Detection of Monoamine  
Neurotransmitters in Live Cells**

Kallol Bera<sup>§</sup>, Anand Kant Das<sup>§</sup>, Ananya Rakshit<sup>§</sup>, Bidyut Sarkar<sup>§</sup>,

Anoop Rawat<sup>§</sup>, Barun Kumar Maity, Sudipta Maiti<sup>§,\*</sup>


<sup>§</sup>Department of Chemical Sciences, Tata Institute of Fundamental Research, Homi Bhabha Road, Colaba, Mumbai

# Contents

1. Reaction Scheme : Dopamine and OPA adduct.....	3
2. Mass spectra of the dopamine OPA adduct.....	3
2.1 MALDI spectra of the dopamin OPA adduct.....	4
2.1.1. Reaction mechanism for OPA dopamine adducts for product I and II...	5
2.1.2. Plausible reaction mechanism of higher molecular weight adduct formation.....	7
2.2 ESI-MS spectra of the Dopamin OPA adduct.....	8
3. OPA and serotonin reaction.....	9
4.1 MALDI spectra of the serotonin (5-HT) OPA adduct.....	10
4.2 ESI-MS spectra of the serotonin (5-HT) OPA adduct.....	11
5. Fluorescence spectroscopy of dopamin-OPA and 5HT-OPA adducts.....	13
6. MN9D, RN46A and HEK293T cell culture.....	13
7. Single photon imaging studies of cells.....	13
8. Colocalization experiments (Multiphoton imaging followed by single photon).....	13
9. Mass spectrometry of OPA treated RN46A Cells.....	14

## 1. Reaction scheme :Dopamine (DA) and OPA adduct

Equivalent concentration of dopamine and OPA are mixed in 10 mM Thomson's buffer (pH=7.4) at room temperature (25 °C) for 30 min in stirring condition.




**Figure S1:** A plausible adducts formation scheme upon reaction with dopamine (DA) and OPA.

## 2. Mass spectra of the dopamine (DA) OPA adduct

1 mM of OPA (in 500  $\mu\text{L}$  10 mM Thomson's buffer, pH = 7.4) and 1 mM of DA: in 500  $\mu\text{L}$  Thomson's buffer, pH = 7.4) are mixed for 30 minutes then the reaction mixture was desalted using ziptip and the reaction products were dissolved in acetonitrile:water (1:1). This mixture was used for confirming possible products using MALDI-TOF mass spectrometer (Bruker model ultrafleXtreme) and ESI mass spectroscopic technique.

## 2.1 MALDI spectra of the dopamine (DA) OPA adduct



**Figure S2:** MALDI TOF mass spectra of dopamine-OPA (DA-OPA) adducts.


$m/z = 250.946$  indicate formation of I ( $m/z_{\text{cal}} = 251.09$ ),

$m/z = 388.328$  indicate formation of II ( $m/z_{\text{cal}} = 367.00$ )

We didn't get the signature of the product III in the MALDI spectrum. However, we have found higher molecular mass which is possibly due to polymerization of DA-OPA adducts through catechol moiety. These polymerization reactions are described in the next section.

### 2.1.1. Reaction mechanism for OPA dopamine adducts for product I and II


#### Product I:



**Figure S3:** Possible reaction mechanism for the product I

The formation of this product I is already confirmed by MALDI (Figure S2)


**Product II:**




**Figure S4:** Possible reaction mechanism for the product II

The formation of this product II is confirmed by MALDI (Figure S2)

### 2.1.2. Plausible reaction mechanism of higher molecular weight adduct formation



**Figure S5:** Possible reaction mechanism for higher molecular mass



Dimer (I''')

Trimer (I''')

**Figure S6:** Structure for the dimer and trimer.

$m/z = 550.669$  indicate formation of DA-OPA dimer I'' ( $m/z_{cal} = 550.17$ )


$m/z = 619.935$  indicate formation of DA-OPA dimer with  $Na^+$  &  $K^+$  salts I''' ( $m/z_{cal} = 619.08$ )

$m/z = 871.392$  indicate formation of DA-OPA trimer I'''' ( $m/z_{cal} = 871.16$ ).

Such polymerization also support formation of I and II.

Note: In the cellular environment, the above polymerization step could not take place as it is reducing in nature.

## 2.2 ESI-MS spectra of the Dopamine (DA) OPA adduct




**Figure S7:** ESI Mass spectra of dopamine-OPA (DA-OPA) adduct.

$m/z = 249.0$  indicate the formation of I' (oxidized form of I)



### 3. OPA-serotonin reaction:




**Figure S8:** A plausible adducts formation scheme upon reaction with serotonin and OPA. Ox stands for oxidation.

**Note:** Serotonin and OPA can form adducts by similar kind of reaction mechanism.

#### 4.1 MALDI spectra of the serotonin (5-HT) OPA adduct

1 mM of OPA (in 500  $\mu$ L 10 mM Thomson's buffer, pH = 7.4) and 1 mM of 5-HT (in 500  $\mu$ L 10 mM Thomson's buffer, pH = 7.4) are mixed for 30 minutes then the mass spectra are checked of the final products after ziptip(C18) using MALDI-TOF mass spectrometer (Bruker model ultrafleXtreme).



**Figure S9:** MALDI TOF mass spectra of serotonin-OPA (5HT-OPA) adducts.


$m/z = 272.006$  indicate formation of IV' ( $m/z_{cal} = 272.09$ ).

$m/z = 274.034$  indicate formation of IV ( $m/z_{cal} = 274.11$ ).

$m/z = 276.030$  indicate formation of VI ( $m/z_{cal} = 276.13$ ).

$m/z = 293.106$  indicate formation of V ( $m/z_{cal} = 292.12$ ).


## 4.2 ESI-MS spectra of the serotonin (5-HT) OPA adduct



**Figure S10:** ESI Mass spectra of serotonin-OPA (5HT-OPA) adduct.

$m/z = 276$  indicate formation of VI ( $m/z_{\text{cal}} = 276.13$ )

$m/z = 293$  indicate formation of V ( $m/z_{\text{cal}} = 292.12$ )



**Figure S11:** Probable serotonin-OPA (5HT-OPA) adducts and their calculated mass.

**5. Fluorescence spectroscopy of DA-OPA and 5HT-OPA adduct:** Fluorescence spectra are recorded using a Fluoromax equipped with a xenon lamp. Samples (DA-OPA and 5HT-OPA adduct) for emission measurements are contained in quartz cuvettes (1 mL volume). Excitation is provided at 450 nm with excitation and emission slit widths at 5/5 nm. All spectroscopic measurements are performed under two different conditions (in 10 mM PBS buffer, pH 7.4 & pH 5.5).

**6. MN9D, RN46A and HEK cell culture.** MN9D dopaminergic cells and HEK cells of low passage numbers were cultured in DMEM media (Gibco, USA) supplemented with 10% FBS, 50 units/ml penicillin and 50 µg/mL streptomycin (Gibco, USA) at 37 °C under humidified air containing 5% CO<sub>2</sub> in T-25 canted-neck flasks. For RN46A serotonergic cells, DMEM-F12 (1:1) (Gibco, USA) media was used instead of regular DMEM. Media was changed every 48 hours.


**7. Single photon imaging of cells.** For cell imaging studies, all cell types were grown on a poly-L-lysine-coated (0.1 mg/mL) homemade glass coverslip-bottomed Petri dishes. Cell imaging was performed using a confocal microscope (LSM 710, Carl Zeiss, Germany) with a 40X water immersion objective. Fluorescence images were recorded at an interval of 0.5 µm in the longitudinal (Z) direction. Different cell types were incubated with 100 µM of OPA in Thomson's buffer (TB, 146 mM NaCl, 5.4 mM KCl, 1.8 mM CaCl<sub>2</sub>, 0.8 mM MgSO<sub>4</sub>, 0.4 mM KH<sub>2</sub>PO<sub>4</sub>, 0.3 mM Na<sub>2</sub>HPO<sub>4</sub>, 5 mM dextrose, and 20 mM Na-HEPES, pH adjusted to 7.4) for 30 min. Subsequent washes were given with TB to remove the remaining OPA. The cells were then imaged (excitation wavelength, 488 nm). All the salt buffers were purchased from sigma.

**8. Colocalization experiments (multiphoton imaging followed by single photon imaging).** The RN46A serotonergic cells (without OPA treatment) were first imaged using multiphoton excitation to capture the intrinsic fluorescence of serotonin, followed by single photon excitation of OPA-serotonin adduct after OPA treatment and probed for co-localization. Briefly, on a custom modified confocal microscope (LSM 710, Zeiss, Germany) an external self mode-locked Ti:Sapphire femtosecond laser (Mai Tai, Spectra-Physics, CA, USA) is used for multiphoton excitation. Briefly, 740 nm light (from the tune-able range of 690-1040 nm of Mai Tai) were routed to the scan box of the confocal microscope. The excitation laser (mildly diverging) was collimated with a 100 cm convex lens. A combination of a half wave-plate and a polarizer was placed to control the excitation power. The custom modification allowed us to attach our detectors and detector housing to the non-descanned port of the LSM 710 confocal microscope. The excitation was routed to the sample through a dichroic (MBS 690+, Zeiss, Germany) in the confocal scan box to another dichroic (MBS 690+, Zeiss, Germany or 720 nm short pass dichroic, Semrock, USA) placed right below the objective lens (40× water immersion Zeiss, Germany) at 45° angle. The fluorescence collected back by the same objective was separated from the excitation by the dichroic below the objective and sent to the non-descanned port of the microscope. The non-descanned port was custom modified to accommodate the external PMT and process the signal recorded by it. We used a near-UV sensitive PMT (P30A-01, Electron tubes, UK), directly coupled to the side exit port, for better detection of UV fluorophores. The

fluorescence was filtered through a 1 cm thick  $\text{CuSO}_4$  solution to get rid of the residual excitation light. The other dichroic in the scan box is just for guiding the excitation light to the sample. For the single photon excitation, the protocol was followed as mentioned in the previous section.

## 9. Mass Spectrometry with OPA treated RN46A Cells

High-passage RN46A cells are cultured in DMEM-F12 (1:1) media supplemented with 10% FBS, 50 units/ml penicillin, and 50  $\mu\text{g}/\text{mL}$  streptomycin at  $37^\circ\text{C}$  under humidified air containing 5%  $\text{CO}_2$  in T-25 canted-neck flasks. Then the media is washed with Thomson's buffer and incubated with 100  $\mu\text{M}$  OPA (2 ml) for 1 hour. Again the cells are washed thoroughly with Thomson's buffer, sonicated for 1 h to crush the cells and the soup is collected by filtration for mass spectrometry. For control experiment we have followed same experiment without OPA.



**Figure S12:** MALDI-TOF mass spectra of (a) RN46A cell + 100  $\mu\text{M}$  OPA extract and (b) only RN46A cell extract.