

Improved Synthesis of 4-Cyanotryptophan and Other Tryptophan Analogues in Aqueous Solvent Using Variants of TrpB from *Thermotoga maritima*

Christina E. Boville, David K. Romney, Patrick J. Almhjell, Michaela Sieben, Frances H. Arnold*

Division of Chemistry and Chemical Engineering 210-41, California Institute of Technology, 1200 East California Boulevard, Pasadena, California 91125, United States

*E-mail: frances@cheme.caltech.edu, Phone: (626) 395-4162, Fax: (626) 568-8743

Table of Contents

Results of site-saturation mutagenesis libraries	S2
Figure S1.....	S2
Figure S2.....	S3
Figure S3.....	S4
LCMS calibration curves for Chart 1	S5
Figure S4.....	S5
Table S1. HPLC data for Figure 2.....	S6
Table S2. HPLC data for reactions after 1 hour.....	S6
Table S3. HPLC data for Chart 1 and indole	S7
Table S4. Primers for construction of recombination libraries	S8
Table S5. Primers for construction of site-saturation libraries	S8
4-cyanotryptophan NMR spectrum for Scheme 2.....	S9
References.....	S9

Results of screening site-saturation mutagenesis libraries




Figure S1. Site-saturation mutagenesis of residue 30 in Tm9D8. Following the procedure in the Experimental Section, 80 colonies (4-fold oversampling) were assayed for 4-CN-Trp production at 50 °C. Although many variants showed parent-like activity, the top two variants were found to be parent (G30).




Figure S2. Site-saturation mutagenesis of residue 228 in *Tm9D8*. Following the procedure in the Experimental Section, 80 colonies (4-fold oversampling) were assayed for 4-CN-Trp production at 50 °C. The top two variants were identified as parent (S228).




Figure S3. Site-saturation mutagenesis of residue 184 in Tm9D8. Following the procedure in the Experimental Section, 80 colonies (4-fold oversampling) were assayed for 4-CN-Trp production at 50 °C. I184L was found in two of the top variants, but this mutation was not as beneficial as I184F in the rescreen.

LCMS calibration curves for Chart 1




Figure S4. Using an authentic product standard, mixtures of starting material and product at different ratios (9:1, 3:1, 1:1, 1:3, and 1:9) were used to generate a calibration curve. Each mixture was prepared in duplicate with a final concentration of 1 mM in 1:1 1-M aq. HCl/ CH₃CN. Mixtures were analyzed with LCMS at 254 nm and 280 nm (reference 360 nm, bandwidth 100 nm) and 330 nm (no reference wavelength) and were correlated to the actual ratios by a linear relationship. The correlation for 4-nitrotryptophan was published previously.¹ For tryptophan and 5-chloro-7-iodotryptophan, the HPLC yield was approximated by comparing absorption peaks of product and starting material at 280 nm.

Table S1. HPLC data for Figure 2

	Temperature (°C)					
	37	37	50	50	75	75
<i>Tm2F3</i>						
254 nm:	5.68%	5.76%	17.1%	17.3%	15.1%	18.7%
280 nm:	3.19%	3.00%	10.8%	10.0%	7.44%	9.78%
<i>Tm9D8</i>						
254 nm:	30.5%	28.7%	38.0%	36.2%	18.0%	16.5%
280 nm:	16.6%	15.2%	21.9%	19.5%	11.7%	10.8%
<i>Tm9D8*</i>						
254 nm:	69.2%	66.5%	70.1%	71.1%	35.1%	35.9%
280 nm:	48.1%	47.6%	49.9%	48.7%	19.3%	20.5%

	Temperature (°C)					
	37	37	50	50	75	75
<i>Tm2F3</i>						
	Corrected					
254 nm:	8.11%	8.22%	23.2%	23.5%	20.6%	25.3%
280 nm:	10.8%	10.2%	30.7%	29.0%	22.8%	28.5%
<i>Tm9D8</i>						
	Corrected					
254 nm:	39.2%	37.2%	47.5%	45.5%	24.4%	22.5%
280 nm:	42.1%	39.6%	50.5%	46.9%	32.7%	30.7%
<i>Tm9D8*</i>						
	Corrected					
254 nm:	77.1%	74.8%	77.9%	78.7%	44.4%	45.1%
280 nm:	76.7%	76.4%	77.9%	77.1%	46.7%	48.6%

Table S2. HPLC data for reactions after 1 hour^a

Catalyst	HPLC yield			Corrected			Average	Std. Dev.
	254 nm	254 nm	254 nm	#1	#2	#3		
<i>Tm2F3</i> ^b	0.834%	0.800%	0.840%	1.22%	1.17%	1.22%	1.20%	0.03%
<i>Tm9D8</i> ^b	3.7%	3.7%	4.3%	5.4%	5.3%	6.1%	5.6%	0.4%
<i>Tm9D8</i> ^{*c}	3.9%	3.8%	4.3%	5.6%	5.5%	6.1%	5.7%	0.3%

^aMaximum 1000 turnovers. ^bReactions conducted at 50 °C. ^cReactions conducted at 37 °C.

**Table S3. HPLC data for Chart 1 and indole
4-bromotryptophan:**

Catalyst	HPLC yield at indicated wavelength							
	254 nm	254 nm	Average	Corrected	280 nm	280 nm	Average	Corrected
<i>Tm2F3</i>	58.3%	58.1%	58.2%	66.5%	60.9%	59.7%	60.3%	64.6%
<i>Tm9D8</i>	63.2%	65.7%	64.4%	71.9%	67.1%	68.0%	67.6%	71.1%
<i>Tm9D8*</i>	69.5%	68.6%	69.1%	75.8%	72.7%	72.2%	72.4%	75.3%

4-nitrotryptophan:

Catalyst	HPLC yield at indicated wavelength							
	254 nm	254 nm	Average	Corrected	330 nm	330 nm	Average	Corrected
<i>Tm2F3</i>	3.45%	2.51%	2.98%	2.52%	1.09%	1.14%	1.11%	2.58%
<i>Tm9D8</i>	1.74%	2.70%	2.22%	1.87%	1.10%	1.08%	1.09%	2.53%
<i>Tm9D8*</i>	–	–	–	–	–	–	–	–

5-nitrotryptophan:

Catalyst	HPLC yield at indicated wavelength							
	254 nm	254 nm	Average	Corrected	330 nm	330 nm	Average	Corrected
<i>Tm2F3</i>	16.4%	17.9%	17.2%	20.9%	19.3%	19.4%	19.3%	20.9%
<i>Tm9D8</i>	4.28%	4.97%	4.62%	5.68%	4.79%	4.92%	4.86%	5.3%
<i>Tm9D8*</i>	5.57%	5.80%	5.68%	6.98%	6.77%	6.97%	6.87%	7.5%

5-bromo-7-fluorotryptophan:

Catalyst	HPLC yield at indicated wavelength							
	254 nm	254 nm	Average	Corrected	280 nm	280 nm	Average	Corrected
<i>Tm2F3</i>	40.8%	42.1%	41.4%	34.3%	46.8%	48.1%	47.4%	34.3%
<i>Tm9D8</i>	96.8%	97.1%	96.9%	94.1%	97.1%	97.7%	97.4%	94.4%
<i>Tm9D8*</i>	94.4%	95.4%	94.9%	91.5%	95.1%	96.4%	95.7%	91.7%

5-chloro-7-iodotryptophan:

Catalyst	HPLC yield		
	280 nm	280 nm	Average
<i>Tm2F3</i>	2.3%	2.2%	2.3%
<i>Tm9D8</i>	19.7%	20.6%	20.2%
<i>Tm9D8*</i>	18.9%	17.8%	18.4%

tryptophan:

Catalyst	HPLC yield		
	280 nm	280 nm	Average
<i>Tm2F3</i>	99.3%	97.3%	98.3%
<i>Tm9D8</i>	90.6%	88.3%	89.4%
<i>Tm9D8*</i>	79.0%	79.7%	79.3%


Table S4. Primers for construction of recombination libraries

Fragment	Forward Primer	Reverse Primer
Fragment 1: <i>Nde</i> I to E30	GAAATAATTTTGTTTAACTTTAAGA AGGAGATATACATATG	CTCATCTTTCATGATTYCTTCGTA CGCAGCTTC
Fragment 2: E30 to I184	GAAGCTGCGTACGAAGRAATCAT GAAAGATGAG	ACCAACCACAGAGCCGTWCACGT AATAGGTGGT
Fragment 3: I184 to G228	ACCACCTATTACGTGWTCGGCTC TGTGGTTGGT	AGCGTTAGAACCACCGCYCACGC ACGCAACGAT
Fragment 4: G228 to <i>Xho</i> I	ATCGTTGCGTGCGTGRGCGGTGG TTCTAACGCT	GCCGGATCTCAGTGGTGGTGGTG GTGGTGCTCGAG
Gene Assembly	GAAATAATTTTGTTTAACTTTAAGA AGGAGATATACATATG	GCCGGATCTCAGTGGTGGTGGTG GTGGTGCTCGAG

Table S5. Primers for construction of site-saturation libraries

target <i>Tm</i> 9D8 residue	forward primer	reverse primer
E30	TGCGTACGAAXXXATCATGAAAGATGAG	GCTTCCAGTTCTTCCAGAG
I184	CTATTACGTGXXXGGCTCTGTGGTTGG	GTGGTCTGCAGGTTGGTA
G228	TGCGTGCGTGXXXGGTGGTTCTA	ACGATGTAGTCCGGCAGAC

The codon indicated as XXX is replaced with NDT, VHG, or TGG.

4-cyanotryptophan NMR spectrum for Scheme 2**References**

- (1) Romney, D. K.; Murciano-Calles, J.; Wehrmüller, J. E.; Arnold, F. H. Unlocking Reactivity of TrpB: A General Biocatalytic Platform for Synthesis of Tryptophan Analogues. *J. Am. Chem. Soc.* **2017**, 139 (31), 10769–10776.